




Product Catalogue

Quality, Service and Value





▶ Panther Medical Canada

History

Panther Medical Canada Inc is a welcome addition to the Canadian Healthcare Supplier Landscape and is an example of Canadians working to reshape healthcare by delivering missing quality, service and value.

By partnering with **Panther Healthcare of China**, the 3rd largest global manufacturer open & endoscopic stapler products, Panther Medical Canada will deliver significant value to Canadian Healthcare.

Panther Medical Canada also partnered with **Purple Surgical of the United Kingdom** and now has ability to deliver a wide range of Surgical Instrumentation.

Panther as well as Purple products are CE Mark Approved, FDA Licensed and Health Canada Licensed, are being used in over 65 Countries WorldWide and are now available in Canadian ORs for clinical use and savings.

Formed in July 2018 by two successful Canadian Healthcare Entrepreneurs, Mario Cortis & David Goossens, Panther Medical Canada plans to partner with Canadian Hospitals, Surgeons and all stakeholders to supply high quality open & laparoscopic surgical solutions.

Panther Medical Canada aims to deliver training, quality and service to Canadian Surgeons and OR Teams, while disrupting pricing inefficiencies in the Canadian Surgical Market.



OUR MISSION

Provide innovative and value healthcare products serving surgical communities



OUR VISION

Source high quality products from best manufacturers across the world and deliver **Quality, Service and Value** to Operating Rooms across Canada



DIRECTIVE PRINCIPLES

- Live Quality
- Execute Lean
- Keep it simple
- Help Society

Contents

1... Staplers, Panther Healthcare, China

Endo Staplers

- 2... Endo Linear Cutter
- 2... Endo Articulating Reloads



Open Staplers

- 3... Linear Cutter and Reloads
- 4... Linear Stapler and Reloads
- 5... Circular Stapler
- 5... Hemorrhoidal Stapler



6... Laparoscopic Instruments, Purple Surgical, UK

Access

- 7... Ultimate Auto Shield Trocars - Bladed
- 8... Ultimate Dilating Tip Trocars - Bladeless
- 9... Ultimate Hasson Trocars
- 9... Ultimate Universal Cannula
- 10... Essential Optical Trocars



Laparoscopic Essentials

- 14... 5mm Blunt Dissector
- 14... 5mm Aspiration Needle
- 14... Fog Free
- 15... Ultimate Veress Needle
- 15... Insufflation Filter and Tubing Sets
- 15... Smoke Filter
- 16... Light Handle Cover
- 16... Instrument Pouch



Laparoscopic Instruments

- 11... Essentials Hand Instrument Range
- 11... Essentials Plus Hand Instrument Range
- 12... Ultimate Suction/Irrigation Sets
- 12... Ultimate Suction/Irrigation Shafts
- 12... Dual Chamber Irrigation Pump
- 13... Ultimate Endo Retrieval Pouch
- 13... Secure Catch Tissue Retrieval Pouch
- 13... Tissue Retrieval Pouch



17... Canada Key Contacts

18... Standard Terms and Conditions

Ask Us About Procedure Kits

- Customized kits as per surgeon or hospital specifications
- Available in breather bags, trays and blister packs

Product Categories

Endo Staplers

Open Staplers

Access

Laparoscopic Instruments

Laparoscopic Essentials

Panther Healthcare, China

Panther Healthcare is the 3rd largest manufacturer of surgical staplers and is globally present in more than 70 countries. Panther Healthcare is focused on surgical stapling and minimally invasive surgery products. Our current portfolio of high quality open and endoscopic Surgical Staplers are serviced by a dedicated team of Clinical Product Specialists.

The Panther HealthCare R&D Team has 4 major product focuses including Endo Staplers, Linear and Circular Cutters and Linear Staplers. Coming soon to Canada are Powered Endo Staplers and many other innovations to help Canadian surgeons deliver better outcomes.



Endo Linear Cutter

The Endo Linear Cutter is available in 3 different shaft lengths- Short, Standard and XL that enable surgeons to perform multiple applications and cover all surgery procedures. The Panther Endo Linear Stapler simultaneously cuts and divides tissue while placing two triple staggered rows of titanium staples in the tissue. The universal stapler handle is designed to host all cartridges of all sizes for up to 25 fires per patient.

Endo Liner Cutter



Product Code	Description	Shaft Length
CEAB30	Endo Linear Cutter (Standard Handle)	6cm
CEAB45	Endo Linear Cutter (Short Handle)	16cm
CEAB60	Endo Linear Cutter (Extra Long Handle)	26cm

Articulating Reload Units



Product Code	Staple Length	Staple Height
CADB-30T	30mm	2.5mm
CADB-30N	30mm	3.5mm
CADB-45T	45mm	2.5mm
CADB-45N	45mm	3.5mm
CADB-45S	45mm	4.2mm
CADB-45D	45mm	4.8mm
CADB-45R	45mm	5.0mm
CADB-60T	60mm	2.5mm
CADB-60N	60mm	3.5mm
CADB-60S	60mm	4.2mm
CADB-60D	60mm	4.8mm
CADB-60R	60mm	5.0mm

Curved Tip Articulating Reload Units



Product Code	Staple Length	Staple Height
CADC-30T	30mm	2.5mm
CADC-45T	45mm	2.5mm

Linear Cutter Stapler

The Panther Linear Cutter Stapler places two double staggered rows of titanium staples while simultaneously cutting and dividing the tissue in between. The staplers and reloads are available in 3 different lengths – 60mm, 80mm and 100mm and 2 staple sizes- 3.8mm (Blue) and 4.8mm (Green). The stapler can be reloaded up to 7 times for total of 8 firings in a single procedure.

Linear Cutter Stapler



Product Code	Description	Staple Length	Staple Height
SSAB-60 3.8mm	Linear Cutter Stapler	60mm	3.8mm
SSAB-60 4.8mm	Linear Cutter Stapler	60mm	4.8mm
SSAB-80 3.8mm	Linear Cutter Stapler	80mm	3.8mm
SSAB-80 4.8mm	Linear Cutter Stapler	80mm	4.8mm
SSAB-100 3.8mm	Linear Cutter Stapler	100mm	3.8mm
SSAB-100 4.8mm	Linear Cutter Stapler	100mm	4.8mm

Cartridges for Linear Cutter Stapler



Product Code	Staple Length	Staple Height
SADB-60N	60mm	3.8mm
SADB-60D	60mm	4.8mm
SADB-80N	80mm	3.8mm
SADB-80D	80mm	4.8mm
SADB-100N	100mm	3.8mm
SADB-100D	100mm	4.8mm

Open Linear Stapler

The Panther Linear Stapler places a double staggered row of titanium staples. The unique color indicator safety feature shows different colors for different firing phases thus significantly reducing the possibility of error. The stapler reloads are available in lengths of 30mm, 45mm, 60mm and 90mm with 3 different staple heights- 2.5mm (White), 3.5mm (Blue) and 4.8mm (Green). The stapler can be reloaded up to 7 times for total of 8 firings in a single procedure.

Open Linear Stapler



Product Code	Description	Staple Length	Staple Height
FLSLF30-3.5	Linear Stapler	30mm	3.5mm
FLSLF30-4.8	Linear Stapler	30mm	4.8mm
FLSLF45-3.5	Linear Stapler	45mm	3.5mm
FLSLF45-4.8	Linear Stapler	45mm	4.8mm
FLSLF60-3.5	Linear Stapler	60mm	3.5mm
FLSLF60-4.8	Linear Stapler	60mm	4.8mm
FLSLF90-3.5	Linear Stapler	90mm	3.5mm
FLSLF90-4.8	Linear Stapler	90mm	4.8mm

Cartridges for Linear Cutter Stapler



Product Code	Staple Length	Staple Height
FLSLH3.5A	30mm	3.5mm
FLSLH4.8A	30mm	4.8mm
FLSLH3.5B	45mm	3.5mm
FLSLH4.8B	45mm	4.8mm
FLSLH3.5C	60mm	3.5mm
FLSLH4.8C	60mm	4.8mm
FLSLH3.5E	90mm	3.5mm
FLSLH4.8E	90mm	4.8mm

Circular Stapler

The Panther disposable Circular Stapler places a circular, double staggered row of titanium staples while resecting the tissue and creating a circular anastomosis. The unique combination of variable staple height with tilt-top anvil delivers secure anastomosis over a wide range of tissue, with easy and safe removal post firing. The Circular Staplers are available in 6 different diameters – 21mm, 24mm, 26mm, 29mm, 32mm and 34mm.

Circular Stapler



Product Code	Description	Closed Staple Height Range	Outer Diameter	Staple Height
FCSME24	Circular Stapler with tilt top	1.4-2.2mm	24mm	4.5mm
FCSME26	Circular Stapler with tilt top	1.6-2.4mm	26mm	4.8mm
FCSME29	Circular Stapler with tilt top	1.6-2.4mm	29mm	4.8mm
FCSME32	Circular Stapler with tilt top	1.8-2.6mm	32mm	5.0mm
FCSLWBE21	Circular Stapler	1.4-2.2mm	21mm	4.5mm
FCSLWBE34	Circular Stapler	1.8-2.6mm	34mm	5.0mm

Hemorrhoidal Staplers

The Panther Hemorrhoidal / Prolapse Circular Stapler is designed for rectal prolapse and hemorrhoid disease. The unique transparent staple casing enables surgeons to easily see the amount of tissue drawn in the casing. Each stapler has a patented dilator with light source, which provides constant illumination of the surgical area. Fine details such as the purse string are clearly more visible in fine detail than traditional PPH dilators. The staplers also boast a patented inside safety trigger to prevent misfirings.

Circular Stapler



Product Code	Description	Outer Diameter
FCSSWAE32	Disposable Hemorrhoidal Stapler	32mm
FCSSWAE33	Disposable Hemorrhoidal Stapler	33mm
FCSSWAE34	Disposable Hemorrhoidal Stapler	34mm

Purple Surgical, UK

Purple Surgical is the UK's leading manufacturer of laparoscopic instruments and devices.

Through local sourcing and advanced British manufacturing of surgical products, Purple Surgical has complete control over the design and evolution of our products. Through the partnership with Purple Surgical, Panther Medical Canada is able to adapt quickly to the constantly evolving healthcare market, while offering high quality, single-use, and cost-effective solutions.

Purple Surgical actively invests in both cutting edge technology and equipment, as well as in the people who design, develop and manufacture all of our high quality surgical products.



Ultimate AutoShield Trocars - Bladed

The Purple Surgical Ultimate AutoShield Trocar Range features an ultra-sharp, linear bi-lateral blade for simple, minimal force insertion into the abdominal cavity.

Ultimate AutoShield Trocars - Bladed



Product Code	Description	Cannula	Pack Qty.
PS3575ULT	5mm x 70mm Trocar	Ribbed	5 / Box
PS3565ULT	5mm x 100mm Trocar	Ribbed	5 / Box
PS3560ULT	10mm x 100mm Trocar	Ribbed	5 / Box
PS3562ULT	12mm x 100mm Trocar	Ribbed	5 / Box
PS3576ULT	5mm x 70mm Trocar + Additional Cannula	Ribbed	5 / Box
PS3563ULT	5mm x 100mm Trocar + Additional Cannula	Ribbed	5 / Box
PS3564ULT	10mm x 100mm Trocar + Additional Cannula	Ribbed	5 / Box
PS3561ULT	12mm x 100mm Trocar + Additional Cannula	Ribbed	5 / Box
PS4575ULT	5mm x 70mm Trocar	Smooth	5 / Box
PS4565ULT	5mm x 100mm Trocar	Smooth	5 / Box
PS4560ULT	10mm x 100mm Trocar	Smooth	5 / Box
PS4562ULT	12mm x 100mm Trocar	Smooth	5 / Box
PS4576ULT	5mm x 70mm Trocar + Additional Cannula	Smooth	5 / Box
PS4563ULT	5mm x 100mm Trocar + Additional Cannula	Smooth	5 / Box
PS4564ULT	10mm x 100mm Trocar + Additional Cannula	Smooth	5 / Box
PS4561ULT	12mm x 100mm Trocar + Additional Cannula	Smooth	5 / Box
PS5575ULT	5mm x 70mm Trocar	Threaded	5 / Box
PS5565ULT	5mm x 100mm Trocar	Threaded	5 / Box
PS5560ULT	10mm x 100mm Trocar	Threaded	5 / Box
PS5562ULT	12mm x 100mm Trocar	Threaded	5 / Box
PS5576ULT	5mm x 70mm Trocar + Additional Cannula	Threaded	5 / Box
PS5563ULT	5mm x 100mm Trocar + Additional Cannula	Threaded	5 / Box
PS5564ULT	10mm x 100mm Trocar + Additional Cannula	Threaded	5 / Box
PS5561ULT	12mm x 100mm Trocar + Additional Cannula	Threaded	5 / Box

Ultimate Dilating Tip Trocars - Bladeless

The Purple Surgical Ultimate Dilating Tip Trocar Range separates rather than lacerates the tissue.

Ultimate Dilating Tip Trocars - Bladeless



Product Code	Description	Cannula	Pack Qty.
PS3475ULT	5mm x 70mm Trocar	Ribbed	5 / Box
PS3465ULT	5mm x 100mm Trocar	Ribbed	5 / Box
PS3460ULT	10mm x 100mm Trocar	Ribbed	5 / Box
PS3462ULT	12mm x 100mm Trocar	Ribbed	5 / Box
PS3476ULT	5mm x 70mm Trocar + Additional Cannula	Ribbed	5 / Box
PS3463ULT	5mm x 100mm Trocar + Additional Cannula	Ribbed	5 / Box
PS3464ULT	10mm x 100mm Trocar + Additional Cannula	Ribbed	5 / Box
PS3461ULT	12mm x 100mm Trocar + Additional Cannula	Ribbed	5 / Box
PS4475ULT	5mm x 70mm Trocar	Smooth	5 / Box
PS4465ULT	5mm x 100mm Trocar	Smooth	5 / Box
PS4460ULT	10mm x 100mm Trocar	Smooth	5 / Box
PS4462ULT	12mm x 100mm Trocar	Smooth	5 / Box
PS4476ULT	5mm x 70mm Trocar + Additional Cannula	Smooth	5 / Box
PS4463ULT	5mm x 100mm Trocar + Additional Cannula	Smooth	5 / Box
PS4464ULT	10mm x 100mm Trocar + Additional Cannula	Smooth	5 / Box
PS4461ULT	12mm x 100mm Trocar + Additional Cannula	Smooth	5 / Box
PS5475ULT	5mm x 70mm Trocar	Threaded	5 / Box
PS5465ULT	5mm x 100mm Trocar	Threaded	5 / Box
PS5460ULT	10mm x 100mm Trocar	Threaded	5 / Box
PS5462ULT	12mm x 100mm Trocar	Threaded	5 / Box
PS5476ULT	5mm x 70mm Trocar + Additional Cannula	Threaded	5 / Box
PS5463ULT	5mm x 100mm Trocar + Additional Cannula	Threaded	5 / Box
PS5464ULT	10mm x 100mm Trocar + Additional Cannula	Threaded	5 / Box
PS5461ULT	12mm x 100mm Trocar + Additional Cannula	Threaded	5 / Box

Ultimate Hasson Trocars

The Purple Surgical Ultimate Hasson Trocar range is designed for surgeons who prefer the security of the open, cut down technique for the primary port.

Ultimate Hasson Trocars



Product Code	Description	Cannula	Pack Qty.
PS3260ULT	10mm x 100mm Trocar	Ribbed	5 / Box
PS3262ULT	12mm x 100mm Trocar	Ribbed	5 / Box
PS4260ULT	10mm x 100mm Trocar	Smooth	5 / Box
PS4262ULT	12mm x 100mm Trocar	Smooth	5 / Box

Ultimate Universal Cannula

The Purple Surgical Ultimate Universal Cannula is compatible with the Ultimate AutoShield, Ultimate Dilating Tip and Ultimate Hasson Trocars.

Ultimate Universal Cannula



Product Code	Description	Cannula	Pack Qty.
PS3577ULT	5mm x 70mm Cannula	Ribbed	5 / Box
PS3567ULT	5mm x 100mm Cannula	Ribbed	5 / Box
PS3568ULT	10mm x 100mm Cannula	Ribbed	5 / Box
PS3566ULT	12mm x 100mm Cannula	Ribbed	5 / Box
PS4577ULT	5mm x 70mm Cannula	Smooth	5 / Box
PS4567ULT	5mm x 100mm Cannula	Smooth	5 / Box
PS4568ULT	10mm x 100mm Cannula	Smooth	5 / Box
PS4566ULT	12mm x 100mm Cannula	Smooth	5 / Box
PS5577ULT	5mm x 70mm Cannula	Threaded	5 / Box
PS5567ULT	5mm x 100mm Cannula	Threaded	5 / Box
PS5568ULT	10mm x 100mm Cannula	Threaded	5 / Box
PS5566ULT	12mm x 100mm Cannula	Threaded	5 / Box

Essential Optical Trocars

See things 360° clearer with enhanced optically guided trocar entry. The Purple Surgical Essentials Optical Trocar range allows for the safe insertion under direct vision eliminating the need for blind entry access.

Essential Optical Trocars



Product Code	Description	Cannula	Pack Qty.
PS5525	5mm x 100mm Trocar	Threaded	5 / Box
PS5526	12mm x 100mm Trocar	Threaded	4 / Box
PS5527	5mm x 100mm Cannula	Threaded	5 / Box
PS5528	12mm x 100mm Cannula	Threaded	4 / Box

Laparoscopic Essentials Instrument Range

The Purple Surgical Laparoscopic Essentials Instrument range is a reliable, highly cost competitive range of monopolar and bipolar 5mm instruments that offer both value and quality in a single-use device.

Laparoscopic Essentials Instrument Range



Product Code	Description	Colour	Pack Qty.
PS4551	Metzenbaum Monopolar Scissors 5mm x 33cm without Ratchet	Purple	5 / Box
PS4581	Mini-Metzenbaum Monopolar Scissors 5mm x 33cm without Ratchet	Purple	5 / Box
PS4552	Maryland Dissecting Forceps 5mm x 33cm without Ratchet	Blue	5 / Box
PS4553	Straight Grasping Forceps 5mm x 33cm without Ratchet	Green	5 / Box
PS4556	Straight Grasping Forceps 5mm x 33cm with Ratchet	Green	5 / Box
PS4555	Lap Clinch 5mm x 33cm without Ratchet	Green	5 / Box
PS4554	Lap Clinch 5mm x 33cm with Ratchet	Green	5 / Box
PS4558	Rat Tooth Lap Clinch 5mm x 33cm with Ratchet	Green	5 / Box
PS4587	Intestinal Forceps 5mm x 33cm without Ratchet	Black	5 / Box
PS4557	Intestinal Forceps 5mm x 33cm with Ratchet	Black	5 / Box
PS4559	Babcock 5mm x 33cm with Ratchet	Black	5 / Box
PS4582	Bipolar Maryland Dissecting Forceps 5mm x 33cm without Ratchet	Yellow	5 / Box

Laparoscopic Essentials Plus Instrument Range

The Purple Surgical Laparoscopic Essentials Plus range is single-use laparoscopic instruments with an on/off ratchet functionality for ultimate surgeon comfort and control.

Laparoscopic Essentials Plus Instrument Range



Product Code	Description	Colour	Pack Qty.
PS4551P	Metzenbaum Monopolar Scissors 5mm x 33cm without Ratchet	Purple	5 / Box
PS4552P	Maryland Dissecting Forceps 5mm x 33cm with on/off Ratchet	Blue	5 / Box
PS4554P	Lap Clinch 5mm x 33cm with on/off Ratchet	Green	5 / Box
PS4558P	Rat Tooth Lap Clinch 5mm x 33cm with on/off Ratchet	Green	5 / Box
PS4556P	Straight Grasping Forceps 5mm x 33cm with on/off Ratchet	Green	5 / Box
PS4557P	Intestinal Forceps 5mm x 33cm with on/off Ratchet	Black	5 / Box
PS4559P	Babcock 5mm x 33cm with on/off Ratchet	Black	5 / Box

Ultimate Suction / Irrigation Sets

The Purple Surgical Ultimate Laparoscopic Suction / Irrigation System has been designed in consultation with leading surgeons. The ergonomically designed device ensures a constant supply of pressurised irrigation fluid and incorporates a leak free seal. The unique 360 degree rotating piston head offers better maneuverability.

Ultimate Suction / Irrigation Sets



Product Code	Description	Pack Qty.
PS3850ULT	Dual Spike 5mm, 32cm Shaft with Tubing, Single-use	10 / box
PS3851ULT	Single Spike 5mm, 32cm Shaft with Tubing, Single-use	10 / box

Ultimate Suction / Irrigation Shafts

The Purple Surgical Ultimate Suction/Irrigation Shafts are compatible with the Ultimate Suction/Irrigation System.

Ultimate Suction / Irrigation Shafts



Product Code	Description	Pack Qty.
PS3852ULT	10mm, Standard Length (32cm) Shaft, Reusable	1 / box
PS3853ULT	10mm, Standard Length (32cm) Shaft, Single-use	10 / box
PS3854ULT	5mm, Extended Length (44cm) Shaft, Single-use	10 / box

Dual Chamber Irrigation Pump

The Purple Surgical Dual Chamber Irrigation Pump is a reusable pump that provides a controlled, uninterrupted high flow of irrigation fluid with adjustable pressure.

Dual Chamber Irrigation Pump



Product Code	Description	Pack Qty.
PS3871ULT	Dual Chamber 2 x 1 Litre	1 / box

Ultimate Endo Retrieval Pouch

The Purple Surgical Ultimate Endo Retrieval Pouch is designed for the easy capture and removal of resected tissue during laparoscopic surgery.

Ultimate Endo Retrieval Pouch



Product Code	Description	Pack Qty.
PS3442ULT	10mm Introducer with 90mm x 145mm, 250ml Pouch	5 / box
PS3443ULT	10mm Introducer with 130mm x 180mm, 580ml Pouch	5 / box
PS3444ULT	15mm Introducer with 190mm x 230mm, 1150ml Pouch	5 / box

Secure Catch Tissue Retrieval Pouch

The Purple Surgical Secure Catch Tissue Retrieval Pouches are designed for the quick, simple, secure capture and removal of resected tissue such as ovaries, gall bladder, fibroids, myomas, omentum, spleen, appendices and bowel.

Secure Catch Tissue Retrieval Pouch



Product Code	Description	Pack Qty.
PS3402	Without Tail - 75mm x 150mm, 250ml Pouch	25 / box
PS3403	Without Tail - 125mm x 175mm, 550ml Pouch	25 / box
PS3404	Without Tail - 175mm x 225mm, 1600ml Pouch	25 / box
PS3405	With Tail - 75mm x 150mm, 250ml Pouch	25 / box
PS3406	With Tail - 125mm x 175mm, 550ml Pouch	25 / box

Tissue Retrieval Pouch

The Purple Surgical Tissue Retrieval Pouches offer the most cost effective solution for secure tissue capture.

Tissue Retrieval Pouch



Product Code	Description	Pack Qty.
PS3400	10mm Introducer with 75mm x 150mm, 120ml Pouch	10 / box
PS3401	10mm Introducer with 125mm x 175mm, 500ml Pouch	10 / box

5mm Blunt Dissector

The Purple Surgical 5mm Blunt Tip Dissector allows for the safe, atraumatic manipulation and dissection of tissue during abdominal investigation.

5mm Blunt
Dissector



Product Code	Description	Pack Qty.
PS3885	5mm Blunt Dissector - Triple Pack	3x10 / box
PS3886	5mm Blunt Dissector - Single Pack	25 / box

5mm Aspiration Needle

The Purple Surgical 5mm Aspiration Needle is designed for the simple, safe aspiration of fluids from cysts and organs.

5mm Aspiration
Needle



Product Code	Description	Pack Qty.
PS3895	5mm Aspiration Needle	10 / box

Fog Free

The Purple Surgical Fog Free Solution prevents the formation of condensation on the distal lens of the laparoscope.

Fog Free



Product Code	Description	Pack Qty.
PS3500	6ml Solution with Applicator Sponge	25 / box

Ultimate Veress Needle

The Purple Surgical Ultimate Veress Needles are designed for the simple, safe, effective creation of the pneumoperitoneum.

Ultimate Veress Needle



Product Code	Description	Pack Qty.
PS3710ULT	120mm Length	25 / box
PS3760ULT	150mm Length	25 / box

Insufflation Filter and Tubing Sets

The Purple Surgical Insufflation Filter & Tubing Set effectively prevents cross contamination between patient and insufflator.

Insufflation Filter and Tubing Sets



Product Code	Description	Pack Qty.
PS3600	With Luer Lok™ Fitting	25 / box
PS3601	With O Ring Fitting (Smith & Nephew)	25 / box
PS3602	With Push Fit Fitting (Stryker)	25 / box
PS3603	With 15mm Fitting (Wolf / Storz)	25 / box
PS3605	With 22mm Fitting (Linvatec)	25 / box

Smoke Filter

The Purple Surgical Smoke Filter enables safe, simple and effective removal of surgical smoke. Diathermy and laser technologies extensively used in laparoscopic surgery result in smoke plume generation. Surgical smoke not only inhibits the surgeons' view of the operative site but also, potentially contains biological pathogens (bacteria, viruses and cancerous cells) as well as chemical toxins.

Smoke Filter



Product Code	Description	Pack Qty.
PS3604	Smoke Filter with Thumbwheel Clamp	25 / box

Light Handle Cover

The Purple Surgical Light Handle Covers are single use, sterile cover that fit easily and securely over our universal operating theatre light handle, which is supplied on a 'free loan' basis.

Light Handle
Cover



Product Code	Description	Pack Qty.
PS4611	Light Handle Cover - Single Pack	100 / box

Instrument Pouch

The Purple Surgical Instrument Pouch enables the secure, convenient housing of instruments close to the surgical site, thus minimising the necessity for hand-to-hand transfer.

Instrument Pouch



Product Code	Description	Pack Qty.
PS4050	Instrument Pouch	25 / box

Canada - Key Contacts Panther Medical Canada

Territory	Name	Phone	Email
British Columbia	Tim Matthews	250-818-7755	tmathews@panthermedical.ca
British Columbia	Rob Horne	604-787-5185	rhorne@panthermedical.ca
Alberta, Saskatchewan, Manitoba	Elvis Millevoi	780-906-7671	emillevoi@panthermedical.ca
SW Ontario	Sean Robertson	519-572-4613	seanrobertson@panthermedical.ca
Quebec and Atlantic Canada	Mouhanad Khouhani	514-446-4465	mkhoulani@panthermedical.ca
National Sales Director	Amandeep Gambhir	416-829-0277	agambhir@panthermedical.ca
National Clinical Director	Lisa DeAngelis	416-991-9444	ldeangelis@panthermedical.ca
General Manager	Mario Cortis	416-274-4705	agambhir@panthermedical.ca

Email: info@panthermedical.ca

Toll Free: 1-833-293-8971

Fax: 1-855-329-5302

Standard Terms and Conditions

ORDER PLACEMENT:

All purchase orders must be emailed as a Purchase Order to orders@panthermedical.ca
All purchase orders must be in typed format, sent as a PDF file, on purchasing entity letterhead.

ORDER SHIPPING:

All purchase orders received before 2:00pm EST will be shipped out same day by Standard ground Purolator Courier.
All orders will be charged the standard Purolator charge for the shipment as charged by Purolator at time of shipping.
Next day delivery rush orders will be charged the Expedited Purolator fee as charged by Purolator at time of shipping.

PRODUCT RETURNS:

Please contact orders@panthermedical.ca to request a Panther Medical Canada Return Goods Authorization Number(RGA#). Any product that is shipped to Panther Medical Canada without a PMC RGA # issued will not be accepted. When an RGA # is obtained, please ship the product pre-paid to Panther Medical Canada with the RGA # clearly marked on outer package/s containing the product return. All product being returned with an RGA # must be in sterile, re-sellable condition and a copy of the invoice must be attached. 15% restocking fee will be charged if the product is returned within 30 days of original purchase date. No credit will be given if the product is returned after 30 days of the original purchase date.

PAYMENT TERMS:

30 days Net on Invoice from date of invoice. No statement will be issued.
1% discount given if the invoice is paid within 10 days of Invoice issuance date.
1.5% monthly interest on invoices not paid within 30 days of original invoice date.

NOTES:

This policy is subject to change without prior notice.



Quality, Service and Value

Head Office:

785 Westney Rd. S Unit 23

Ajax, ON L1S 7G1

Tel: 1-833-293-8971

Fax: 7-855-329-5302

Email: ask@panthermedical.ca

www.panthermedical.ca

