












Panther Medical Canada Inc. Unique Product Features

Endo Linear Cutter	Product	Description	Feature	Johnson & Johnson™	Medtronic™/Covidien™
	CEAB	Universal Endo Handle	Non stop articulating lever		
	CADA and CADB (and CADC for Thoracic surgery)	Straight and Articulated cartridges for Endo linear cutter stapler	Largest, most complete portfolio of endo reloads in the market: Range of RU staple lengths include: WHITE (2.5mm) up to and including GOLD (4.2mm - as with J&J™) and BLACK (5.0mm - as with Medtronic™/Covidien™). All RU's = SLIP RESISTANT. RU'S & ALL HANDLES FULLY INTERCHANGABLE	Has 4.2mm GOLD but not 5.0mm BLACK	Has 5.0mm BLACK but no 4.2mm GOLD
			Each Replacement Unit (RU) has a NEW stainless blade.	Not available	Available
Circular stapler	Product	Description	Feature	Johnson & Johnson™	Medtronic™/Covidien™
	FCSME	 Circular Stapler	Full range of lumen size circular staplers, 21mm up to 34mm. ALL sizes come with dial-in compression setting (as with J&J™) Unique feature: 29mm with tilt top	Has no tilt top	Has no dial-in compression
Linear Stapler with Knife	Product	Description	Feature	Johnson & Johnson™	Medtronic™/Covidien™
	FLSLF	Linear Stapler with Color indicators All sizes (3.5mm to 4.8mm staple leg length) and all lengths (30mm up to 90mm). particularly desired in Thoracic Surgery) (90mm)	SAFETY COLOR INDICATORS - different colors for different firing phases and reduces the possibility of misoperation. 	90 mm not available	90 mm not available
Linear Cutter	Product	Description	Feature	Johnson & Johnson™	Medtronic™/Covidien™
	SAAB and SADB	Linear Cutter and Replacement Unit	Reload Units have new stainless blade for each cutting	Not available	Available
			Panther's Linear Cutter safety interlock physically prevents users from clamping. Surgical operators will notice the "white flag" and realize they need to reload a fresh unit and a used cartridge will never make it back to the patient.	Ethicon™ has a safety lock out to prevent firing on a used cartridge. But, Surgeons can position the Ethicon™ TLC, clamp in lumen and start to fire before they realize they need to reload.	Unknown
Hemorrhoidal stapler	Product	Description	Feature	Johnson & Johnson™	Medtronic™/Covidien™
	FCSSWAE	Hemorrhoidal Stapler	3 sizes: 32, 33 and 34	Only size 33 available	Only size 33 available
			Transparent cartridge	Not available	Not available

Notes:

Ethicon™ is a Trademark of Ethicon Inc. Inc, a subsidiary of Johnson and Johnson Inc.

Medtronic™ is a Trademark of Medtronic Inc.

Disclaimer - All information is subject to changes.